

Negative Pressure Wound Therapy Dressing Kit Breakdown

Kit Contents	Kits				
	Small	Medium	Large	X-Large	Wooding-Scott
Drain	Flat or 10fr	Flat or Channel	Flat or 19fr	19 fr	Irrigation aspiration
Non-adherent gauze	1 – 3 x 3”	1 – 3 x 8”	3 – 3 x 8”	3 – 3 x 8”	3 – 3 x 8”
Antimicrobial Gauze	1 – package, 4 x 4s	1 – package, 6 x 6s	1 – Gauze roll	2 – Gauze roll	2 – Gauze roll
Transparent film	2 – 4 x 4.75” sheets	2 – 6 x 8” sheets	2 – 8 x 12” sheets	2 – 8 x 12” sheets	2 – 8 x 12” sheets
Connective tubing	1 each	1 each	1 each	1 each	1 each
Saline bullet	1 bullet	2 bullet	3 bullet	4 bullet	4 bullet
Skin prep	1 - packet	1 - packet	2 - packets	2 - packets	3 - packets
Strip paste	1 strip	1 strip	1 strip	1 strip	1 strip
Wound ruler	1 each	1 each	1 each	1 each	1 each
Waterproof tape	2 strips	2 strips	4 strips	4 strips	4 strips